



Blood Glucose Monitoring System Gold-Accu

Specification

Blood sample: Fresh capillary whole blood, venous whole blood, arterial whole blood

Calibration sample: Venous plasma

Principle: FAD-GDH

Test range: 1.1-33.3 mmol/L(20-600mg/dL)

Test time: 5 seconds

Blood volume: Less than 0.8 ul

HCT: 20-70%

Expiry date: 24 month

Memory: 500sets blood glucose test values,
100 sets control solution test values

Battery life: 1000 tests

Discard the remaining strip 6 month after the first opening date



Product advantage

- High accuracy, Long validity
- No coding, Easy to use
- Wide HCT range, Suitable for hospital



Sinocare Inc.

Add: No.2F5, Guyuan Road, Hi-Tech Zone, Changsha, Hunan, China
Tel: 86-731-89935581 89935582 Fax: 86-731-89825189
Email: info@sinocare.com
www.sinocare.com



SINOCARE