

PRODUCT PROFILE

H-100 Urine Analyzer

Product Characteristics:

- Adopting the advanced super-high luminosity cold light source with 4 -wavelength, which improves the sensitivity, accuracy, specificity, and reduces the interference from ambient light
- Adopting automatic waste handling equipment, which avoids cross-contamination between samples
- Automatically rectify the test results influenced by non-specificity, pH, specific gravity, and color
- Built-in thermal printer with high speed and low noise; External stylus printer
- Connectable with urine sediment analyzer
- Users can set an abnormal value flag by themselves
- With Chinese and English options
- International, Regular and Symbol system units available for option



Technical Specification:

- ⌘ Test items: urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocytes, glucose, specific gravity, pH and VC (or micro-albumin)
- ⌘ Test wavelength: 525nm, 572nm, 610nm, 660nm
- ⌘ Test principle: Photoelectric colorimetry
- ⌘ Suitable strips: URITIK H8, H10, H11 and H11-MA urinalysis strips
- ⌘ Test throughput: 120 strips/h or 60 strips/h optional
- ⌘ Data memory: 1000 patient results
- ⌘ Computer interface: RS-232 port; parallel printer interface
- ⌘ Display: 240×64 dot-matrix LCD
- ⌘ Language: Chinese and English
- ⌘ Power supply: 100~240VAC, 50Hz/60Hz
- ⌘ Power: 40VA
- ⌘ Dimensions: 376mm×316mm×170mm
- ⌘ Weight: 3.6kg
- ⌘ Printer: Built-in thermal printer

DIRUI INDUSTRIAL CO., LTD.

Address: 95 Yunhe Street, New & High Tech. Development Zone, Changchun, 130012, P. R. China
Tel: +86 431 85100409 Fax: +86 431 85173354
E-mail: dirui@dirui.com.cn Http://www.dirui.com.cn

DIRUI